Quick Guide

Beyond the PCR technology,
BIOFACT promises the progress for your research.

BioFACT™

2X LAMP Master Mix (DNA&RNA)

Description: BioFACT™ 2X LAMP Master Mix is a Loop-Mediated Isothermal Amplification(LAMP) based solution for faster and simpler detection of viral pathogen, including SARS-CoV-2, influenza, Measles, and other pathogen. Our master mix reagents provide maximum flexibility to optimize and accelerate your viral pathogen research and surveillance.

[Cat. No. RL301-10h, RL301-50h]

Contents	RL301-10h	RL301-50h
BioFACT™ 2X LAMP Master Mix (DNA&RNA)	1.0 mL x 1 ea	1.0 mL x 5 ea

Application:

- Loop-Mediated Isothermal Amplification(LAMP)
- RNA and DNA target detection
- Assay development for pathogen detection

Protocol

1) Add following components for a single 20 µl reaction.

PCR Mixture	Vol. 20 µl	Final Conc.
RNA or DNA Template	X µl	> 10 copies/rxn
BioFACT™ 2X LAMP Master Mix	10 µl	1X
Target-Specific primer Mix	X µl	
Fluorescent Dye(optional)	X µl	
RNase free water	Adjust to final 20 µl	

- 2) After all reagents have been added, mix the reaction completely by pipetting several times.
- 3) If more than one reaction is being run, dispense of the reaction mix for each reaction into PCR tubes or a 96-well PCR plate.
- 4) Add of target RNA or target DNA to the appropriate reaction tubes or wells. Add of RNase free water to the NTC reaction tubes or wells. Mix completely by pipetting.
- 5) Cap tubes or seal plate wells. Centrifuge briefly to collect contents prior to incubation.
- 6) Depending on the detection method chosen, incubate the reaction as follows in a real-time detection instrument, thermocycler, or heat block.
- 7) Incubation the following at 65°C for 30 ~ 60 minutes.
 - For fluorescent assays in real-time detection instruments: Monitor fluorescence using at 1 minutes intervals for 30 ~ 60 minutes.

Important notes:

- ① Use separate workspaces as well as specially dedicated equipment.
- ② Opening of tubes or wells after LAMP in the separate post-reaction area helps eliminate risks of environmental contamination.
- ③ No-template control(NTC) and/or other negative control reactions are strongly recommended, in order to confirm absence of background amplification.
- ④ Closes tubes with NTC before adding target RNA to minimize risks of sample cross-over contamination.

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Please contact us, if you have any question and need help.
T)1670-5695 www.bio-ft.com info@bio-ft.com

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