Quick Guide

Beyond the PCR technology,
BIOFACT promises the progress for your research.

LYOFACT™ LAMP Pre-Mix, (DNA & RNA)

[Cat. No. RL601-096]

Description: LYOFACT™ LAMP Pre-Mix(DNA & RNA) is Loop-Mediated Isothermal Amplification(LAMP) and (RT-LAMP) based Lyophilized kit for faster and simpler detection of viral pathogen, including SARS-CoV-2, influenza, Measles, and other pathogen. Our Pre-Mix reagents provide maximum flexibility to optimize and accelerate your viral pathogen research and surveillance.

| Contents | RL601-096 |
|-----------------------------------|---------------------------|
| LyoFACT™ LAMP Pre-Mix (DNA & RNA) | 96 Tube (8 Strip x 12) |

Application:

- Loop-Mediated Isothermal Amplification(LAMP)
- RNA and DNA target detection
- Assay development for pathogen detection
- Lyophilized Type

Protocol

1) Add following components for a single 20 µl reaction.

| PCR Mixture | Vol. 20 µl | Final Conc. |
|----------------------------------|-----------------------|-----------------|
| DNA /RNA Template | X µl | > 10 copies/rxn |
| LyoFACT™ LAMP Pre-Mix(DNA & RNA) | - µl | 1X |
| Target-Specific Primer Mix | X µl | |
| Fluorescent Dye(optional) | X µl | |
| RNase free water | Adjust to final 20 µl | |

- 2) After all reagents have been added, mix the reaction completely by pipetting several times.
- 3) If more than one reaction is being run, dispense of the reaction mix for each reaction into PCR tubes or a 96-well PCR plate.
- 4) Add of target RNA or target DNA to the appropriate reaction tubes or wells. Add of RNase free water to the NTC reaction tubes or wells. Mix completely by pipetting.
- 5) Cap tubes or seal plate wells. Centrifuge briefly to collect contents prior to incubation.
- 6) Depending on the detection method chosen, incubate the reaction as follows in a real-time detection instrument, thermocycler, or heat block.
- 7) Incubation the following at 65°C for 30 ~ 60 minutes.
 - For fluorescent assays in real-time detection instruments:
Monitor fluorescence using at 1 minutes intervals for 30 ~ 60 minutes.

Important notes:

- ① Use separate workspaces as well as specially dedicated equipment.
- ② Opening of tubes or wells after LAMP or RT-LAMP in the separate post-reaction area helps eliminate risks of environmental contamination.
- ③ No-template control(NTC) and/or other negative control reactions are strongly recommended, in order to confirm absence of background amplification.
- ④ Closes tubes with NTC before adding target DNA/RNA to minimize risks of sample cross-over contamination.

Storage:

- 20±5°C for up to 6 months, or 4±5°C for up to 12 months.
- Store it in an aluminum coated bag and in a dry place.
- Keep the humidity below 65% after opening.

Please contact us, if you have any question and need help.
T)1670-5695 www.bio-ft.com info@bio-ft.com

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