Quick Guide

Beyond the PCR technology, BIOFACT promises the progress for your research.

BioFACT™

Bst DNA Polymerase, Large Fragment

[Cat. No. BP101-16h, BP101-80h]

Contents	BP101-16h	BP101-80h
BioFACT™ Bst DNA Polymerase, Large Fragment (8 unit/μl)	1,600 units	8,000 units
10X Bst Reaction Buffer	0.5 mL x 3 ea	1.0 mL x 3 ea
Each 10 mM dNTP Mix	0.9 mL x 1 ea	0.9 mL x 5 ea

Description: BioFACT™ Bst DNA Polymerase, Large Fragment is the portion of the *Bacillus stearothermophilus* DNA Polymerase protein that contains the 5' → 3' polymerase activity, but lacks 5' → 3' exonuclease activity. Bst DNA polymerase, Large Fragment is prepared from an *E.coli* strain containing the *Bacillus stearothermophilus* DNA polymerase gene, lacking the 5' → 3' exonuclease domain.

Protocol [Typical LAMP Protocol]

1) Add following components for a single 25 μl reaction.

PCR Mixture	Vol. 25 μl	Final Conc.
DNA Template	X μl	> 10 copies/rxn
10X Bst Reaction Buffer	2.5 μl	1X
Each 10 mM dNTP Mix	4.5 μl	each 1.8 mM
FIP/BIP Primer(25X)	1 μl	1.6 μM
Loop F/R Primer(25X)	1 μl	0.4 μM
F3/B3 Primer(25X)	1 μl	0.2 μM
Bst Polymerase, Large Fragment (8 unit/μl)	1 μl	8 unit/rxn
Add D.W to	Adjust to final 25 μl	

2) Incubation the following at 65°C for 30 ~ 60 minutes.

Application :

- Loop-mediated isothermal amplification(LAMP)
- Whole genome amplification(WGA)
- Multiple displacement amplification(MDA)
- Ramification amplification(RAM)
- Random-primed DNA labeling

General Guidelines :

- ① A LAMP Primer Mix can be prepared with all 4 or 6(with Loop) primers. A 25X Primer Mix should contain : 40 μM FIP, 40 μM BIP, 5 μM F3, 5 μM B3, 10 μM LoopF, 10 μM LoopR in TE or water.
- ② Reaction should be setup on ice.
- ③ Running a no-template control is strongly recommended to ensure amplification specificity.
- ④ If optimization is desired, try titrating Bst DNA Polymerase (0.04 ~ 0.32 unit/μl) or changing reaction temperature(50 ~ 68°C)

Quality Control :

- Purity : > 99% on SDS-PAGE
- Endonuclease-free
- Exonuclease-free
- RNase-free

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Please contact us, if you have any question and need help.
T)1670-5695 www.bio-ft.com info@bio-ft.com

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